



Exhaled breath microbiomics and amplicon-based biomonitoring

Dr Sterghios A. Moschos FRSC FIBMS FHEA MRSB

Associate Professor, Cellular and Molecular Sciences

Chief Scientific Officer, PulmoBioMed Ltd.

@docmoschos



Northumbria
University
NEWCASTLE

Overview

- Exhaled Breath Diagnostics.
 - Background.
 - Hypothesis.
- The PBM-xHALE™ approach.
 - Platform.
 - Microbiomics.
- Clinical progress on SARS-CoV-2.



Exhaled Breath Condensate (EBC)

- Breath is 95% hydrated.
 - Volatile Organic Compounds.
 - Vapour & aerosols.
 - Biological molecules.
- Biomarker rich:
 - Lung infections.
 - Liver & metabolic disease.
 - Multiple cancers.



Challenges to clinical use

- Reproducibility.
- Contamination:
 - Saliva.
 - Ambient.
- Sample loss.
- Safety.
- Upper vs deep lung separation.

RTube™



Poor process control

EcoScreen™



Sample lost in black tube
17Kg + weight


Exhaled breath condensates (EBC)

- Boyle's law 1662: $PV = nRT$.
- Dreschel bottle: late 1800's.
- Aerosol physics: 1930's.
- Focus on **isolating the sample**.
- **Eliminate sample loss.**
- **Separate saliva.**




Preliminary data: bacteria detection



Bacteria & fungi die by dry ice condensation



NAATs work with 1-3 breaths (bacteria)




EBC is distinct to saliva



Most similar samples
Come from lung surgery
(60% proteobacteria)

(Sze MA et al. *Am. J. Resp. Crit. Care Med.* 2012)

PBM-HALE™: the platform



EBC collector:

- Volatiles and
- Proteins.
- DNA.
- RNA.
- Lipids.
- Medications.

Solves key problems:

- Reproducibility.
- Contamination.
- Sample loss.
- Safety.

Cold Chain Dependent:

- Uses dry ice powder (CO₂) to collect sample reliably.
- Dry ice replenished every 1 hr from compressed gas cylinder.
- Sample needs on the spot test or frozen transfer to lab.

WO2017153755A1: exhaled breath collector – national phase; WO2019053423A1: cascade impactor array – PCT.

PBM-HALE™: the platform



EBC collector:

- Volatiles and
- Proteins.
- DNA.
- RNA.
- Lipids.
- Medications.

Solves key problems:

- Reproducibility.
- Contamination.
- Sample loss.
- Safety.



Path to removing the cold chain:

- Proprietary coating to remove need for dry ice.
- Stabilisation material to remove freezer storage.

Experiments under way

WO2017153755A1: exhaled breath collector – national phase; WO2019053423A1: cascade impactor array – PCT.


Highly consistent sampling



Sampling linear from 25 sec (screening) to 30 min (discovery)

R^2 range: 0.88 to 0.95, $n = 5$.

No salivary contamination



Saliva enzyme levels at least 5000x above assay limit of detection


n = 5.

No bacterial viability




2 min sampling period
n = 5, blood agar

No bacterial viability



Cooling efficiency lost after ~40 min continuous sampling.
n = 3.

16S Microbiomics




Higher DNA content vs background / kit controls

Detection of *Streptococcus* (FDR q = 0.019) w/out extraction.


Aerosolised viruses

Efficient capture of aerosolized virus; dry ice halves infection risk.
25 infectious virions over 15 min sampling period.



GFP-expressing VSV-pseudotyped lentivirus at MOI 0.01 nebulized using PARI TurboBoy SX and captured using PBM-HALE™. Condensates seeded on 10,000 HEK-293T's and GFP expression measured at 72hrs by FACS, visualized by fluorescent microscopy. Bar = 100 µm.


Consistent [protein]



Requires 5x concentration by lyophilization

n=5

Consistent [RNA]




Increased EBC sampling does not increase [RNA].


EBC volume extracted by Trizol matched with volume obtained after 20 breaths.

n=6

18S RNA




2 step SYBR Gold RT-qPCR (triplicate)
 EBC1-5: RNeasy kit 20 breaths
 EBC6: Trizol 30 min sample
 Nasal = swab.



2 step SYBR Gold RT-qPCR (triplicate)
 -ve 1: No RT control
 -ve 2: no cDNA
 EBC volume normalized to 20 breaths

GAPDH & β -actin RNA



2 step SYBR Gold RT-qPCR (triplicate)

-ve 1: No RT control


-ve 2: no cDNA

EBC volume normalized to 20 breaths

Either low host [RNA] or 18S all fungal.

n=5

Metabolomics



Requires 5x concentration by lyophilization

n=5

Metabolomics

Compound	RMM (g/mol)	RT [min]	Relative ion abundance
1-hexadecyl-glycero-3-phosphate	396.3	1.002	810,094
monoacylglyceride	352.3	1.02	281,866
LysoPA	410.2	1.032	968,316
Palmitoleoylethanolamide	297.3	1.047	187,282
eicosatetraenoate	335.2	1.054	348,544
Linoleamide	279.3	1.061	216,809
Cuscohygrine	224.2	1.067	723,759
N-Decanoylglycine	229.2	1.156	2,612,124
N-Nonanoylglycine	215.2	1.198	1,942,872
cis-3-Hexenyl b-primeveroside	394.2	1.221	160,089
N-Lauroylglycine	257.2	1.923	286,977
N-Undecanoylglycine	243.2	2.072	227,826
phosphatidylethanolamine	837.5	2.388	381,518
Gambogic acid	628.3	2.536	416,778
2-Hexenoylcarnitine	257.2	3.062	994,821
L-argininium	175.1	3.367	502,141
N-Acetylputrescine	130.1	3.519	192,382

Data generated at the Northumbria University Metabolomics Core Service

Compounds detected by MS1

- C6-C24 fatty acids.
- Phospholipids & precursors.
- Glycans.
- Medications.
- Drugs of abuse.
- Dietary compounds.

Additionally:

- 20 multiple HDBM hits.
- 104 novel compounds.

COVID-19 Pilot data.

- <30 min EBC (3-4 ml)
 - n=60 target.
 - n=35 interim.
 - 6-10 centres.
 - In COVID19 red zone/ward.
 - Use local Dx RT-PCR assay.
 - Blinded sample Dx analysis.
- Target population:
 - Any age group/sex.
 - Week 1 of symptoms.
 - Nasal positive (<24 hrs).
 - Focus on symptomatic close contacts of COVID-19 cases.

NO FALSE POSITIVES OR FALSE NEGATIVES TO DATE (n=6 controls).
30-60 sec appears adequate.

SARS-CoV-2 in exhaled aerosols?

- Fast: <1 min sample.
 - Simple: just breathe.
 - >2x safer: process halves infectivity.
 - More reliably: Tuneable sample vs nasal swabs.
-
- For mass screening: by mass production of plastic.
 - Where patients are: using *any* point of need testing system.
 - With current gold standard tests (QuRapID^{®1}, ID NOW[®]).

Next steps

Research

- Eliminate cold chain.
- Detect virus by LFT / aptamers.
- On-board analytical capability.

Development

- Spinout trading 19/10/20.
- MHRA LNO / FDA IDE; EUA.
- Mass production, distribution.
- Adaptations for environmental sampling.

**ALL DEVICES PRODUCED THROUGH CHARITABLE/GOVERNMENT FUNDING
WILL BE MADE FREELY AVAILABLE FOR COVID19 PANDEMIC**

Acknowledgements



Moschos lab

Dr Pep Canyelles-Pericas

Dr John Henderson

Dr Andrew Nelson

Dr William Cheung

Dr Louise Usher

Mr Kavit Shah

Mr Guilherme M.E. Silva

Mr Craig Clements

Mr Declan Gardner

Miss Zoe Hewitson

University of Ulm

Prof. Manfred Frick

Dr. Jan Münch

University of Oxford

Prof. Aris Katzourakis

PulmoBioMed Ltd.

Dr. Pete Hotten (CEO)

Dr. Huw A. Edwards (Chairman)

Mr. Jonathan Brooks (COO).

Mr Saqib Ali



Innovate UK



Consultancy disclosures:



How do I use it?

1. Device use SOP: <https://youtu.be/h6tLt9u-rWU>
2. Lay explanation of use: https://youtu.be/TkQEj-KN_os