

Coronavirus 19

Understanding Breath Aerosol Transmission Risk

Dr Sterghios A. Moschos FRSC FIBMS FHEA MRSB

Chief Scientific Officer, PulmoBiomed Ltd.

Associate Professor, Cellular and Molecular Sciences


@docmoschos



Northumbria
University
NEWCASTLE

Overview

- **Breath for Coronavirus screening.**
 - Scientific evidence.
 - Key unknowns.
- **Exhaled Breath Diagnostics.**
 - Opportunities.
 - Challenge.
 - Competition.
- **The PBM-HALE™ approach.**
 - Platform IP.
 - Supporting key data.
- **The proposition.**



PBM-HALE™ (MVP)

Coronavirus 19 (COVID-19) in breath: Confidence in rationale



Biology

- SARS-CoV-2 binds ACE2 receptor.¹
- ACE2 protein levels highest in lower lung.²
- Aerosols (<5 μm) best to reach lower lung (drug delivery science).³

Pathology

- Disease of the lower lung: respirator need.
- Proposed transmission routes: fomites, droplets (cough, >5 μm), but:
 - Models & data⁴ show transmission without symptoms (no cough!).
 - Aerosol science used in epidemiology out of date.⁵

Coronavirus 19 (COVID-19) in breath: Confidence in rationale




Clinical evidence

- Transmission occurs up to 1 week before symptoms (peak @ -2.9 days).¹
- Virus genome levels max in lower lung samples > nose > throat.²⁻⁵
 - <42% false negative oral swabs.
 - 10-25% false negative nasal swabs.
 - Viable virus levels low in nasal swabs.
 - Nasal detection ~70% days 0-5 from symptoms⁶

Experimental evidence:


- COVID-19 ward aerosol gel traps –ve, but ceiling air vents +ve:⁷ Droplets pulled by gravity, aerosols pulled by air flow.⁸
- Aerosolised virus infectious for 16hrs after mechanical generation.^{9, 10}
- Other coronaviruses can naturally aerosolize (n=3000).¹¹

Coronavirus 19 (COVID-19) in breath: Confidence in rationale



Detection of SARS-CoV-2 in nasal swabs (RNA) or blood (IgX) from symptom onset

Coronavirus 19 (COVID-19) in breath: Confidence in rationale



Our hypothesis

- Disease is a function of amount of virus reaching the lower lung.
- Achieved mainly by breath aerosols (or poor immune system).
- Explains close contact transmission chains.

We need to test breath aerosols for:

- *The amount of virus present (genomes).*
- *Infectivity (viruses).*

Diagnosing from Exhaled Breath Condensates (EBC)



Breath is 95% hydrated:

- Volatile compounds (smells, eg garlic, alcohol).
- Vapour & aerosols.
- Biological molecules.

Health and Disease indicators:

- Lung infections.
- Liver diseases.
- Multiple cancers:
 - Blood.
 - Breast.
 - Brain.

Challenges to clinical use

- Reproducibility.
- Contamination:
 - Saliva.
 - Ambient.
- Sample loss.
- Safety.
- Upper vs deep lung separation.

RTube™



Poor process control

EcoScreen™



Sample lost in black tube
17Kg + weight

PBM-HALE™: the platform



PBM-HALE™
(TLR8 RUO/POC/PON use)

EBC collector:

- Volatiles and
- Proteins.
- DNA.
- RNA.
- Lipids.
- Medications


Solves key problems:

- Reproducibility.
- Contamination.
- Sample loss.
- Safety.

Cold Chain Dependent:

- Uses dry ice powder (CO₂) to collect sample reliably.
- Dry ice replenished every 1 hr from compressed gas cylinder.
- Sample needs on the spot test or frozen transfer to lab.

PBM-HALE™: the platform



PBM-iHALE™
(TLR 5 PON use)

EBC collector:

- Volatiles and
- Proteins.
- DNA.
- RNA.
- Lipids.
- Medications.

Solves key problems:

- Reproducibility.
- Contamination.
- Sample loss.
- Safety.


Path to removing the cold chain:

- Proprietary coating to remove need for dry ice.
- Stabilisation material to remove freezer storage.

Experiments under way


Preliminary data: pathogen DNA

Bacteria & fungi die by dry ice condensation





NB: no colonies on plates, white marks are condensation

COVID-19-like test works with 1-3 breaths for bacteria





EBC is distinct to saliva



Tool device

Closest sample like this comes from lung surgery samples only (60% proteobacteria)
(Sze MA et al. *Am. J. Resp. Crit. Care Med.* 2012)



Prototype: highly consistent sampling



Whether 5 breaths (25 sec; e.g. screening) or 30 min of sampling (e.g. discovery)

R^2 range: 0.88 to 0.95, $n = 5$.

Prototype: no salivary contamination




Saliva enzyme levels below limit of detection in EBC:

At least 5000x less in EBC than in saliva (drool) or device saliva trap levels even after 30 min sampling.

n = 5.


Prototype: No microbial growth due to dry ice condensation (blood agar).



2 min sampling period
(2x target sampling period for COVID-19 screening use).

n = 5

Prototype: microbial growth only after loss of dry ice cooling efficacy




Colonies from dry ice-captured EBC cultured on blood agar.

Lateral contact of sampling tube to dry ice lost ~40 min after continuous sampling.

n = 3.

Prototype: consistent [protein] in EBC




Samples lyophilized and re-constituted in 1/5th of original volume:
No statistically significant difference in concentration over time by micro BCA (data close to LLOD).

No concentration increase anticipated.

n=5

Prototype: consistent [RNA] in EBC




EBC breaths, n=6

EBC 30-100 normalized to 20 breath
sample volume, Trizol extraction


No statistically significant difference in concentration.

260/280 ratios reported per column

Prototype: 18S by PCR in EBC RNA




2 step SYBR Gold RT-qPCR (triplicate)
 EBC1-5: RNeasy kit 20 breaths
 EBC6: Trizol 30 min sample
 Nasal swab.



2 step SYBR Gold RT-qPCR (triplicate)
 -ve 1: No RT control
 -ve 2: no cDNA
 EBC 30-100 normalized to 20 breath sample volume
 P<0.01 Saliva vs EBC, EBC vs -ve, Saliva vs -ve.

Prototype: Human RNA in EBC RNA



2 step SYBR Gold RT-qPCR (triplicate)


-ve 1: No RT control

-ve 2: no cDNA

EBC 30-100 normalized for EBC 20 sample volume

- EBC not classified as human tissue.
 - Human Tissue Act 2004.
 - Cells, DNA, or RNA.
 - EBC explicitly excluded.
- Early data: EBC is 18S+
 - Validation round under way with Taqman® assay.
 - ITS SEQ run planned (human vs fungal).

Prototype: Metabolomics in EBC




Metabolite profile in EBC after 5x concentration by lyophilization, and extraction

n=5

Data generated at the Northumbria University Metabolomics Core Service

John Henderson, William Cheung
U. Northumbria

Prototype: VOCs, C₂₄, lipids, meds, drugs



Compound	RMM (g/mol)	RT [min]	Relative ion abundance
1-hexadecyl-glycero-3-phosphate	396.3	1.002	810,094
monoacylglyceride	352.3	1.02	281,866
LysoPA	410.2	1.032	968,316
Palmitoleylethanolamide	297.3	1.047	187,282
eicosatetraenoate	335.2	1.054	348,544
Linoleamide	279.3	1.061	216,809
Cuscohygrine	224.2	1.067	723,759
N-Decanoylglycine	229.2	1.156	2,612,124
N-Nonanoylglycine	215.2	1.198	1,942,872
cis-3-Hexenyl b-primeveroside	394.2	1.221	160,089
N-Lauroylglycine	257.2	1.923	286,977
N-Undecanoylglycine	243.2	2.072	227,826
phosphatidylethanolamine	837.5	2.388	381,518
Gambogic acid	628.3	2.536	416,778
2-Hexenoylcarnitine	257.2	3.062	994,821
L-argininium	175.1	3.367	502,141
N-Acetylputrescine	130.1	3.519	192,382




Compounds detected by MS1:

- C6-C24 fatty acids.
- Phospholipids & precursors.
- Glycans.
- Medications.
- Drugs of abuse.
- Dietary compounds.

Additionally:

- 20 multiple HDBM hits.
- 104 novel compounds.

Prototype: 16S microbiomics detection of streptococcal mild cough.



Higher DNA content vs background controls.


Detection of Streptococcus (BI FDR $q = 0.019$); amplicon generation w/out extraction.

Detailed work under way to optimize processes, polymerases, pipelines.

John Henderson
Andrew Nelson
Darren Smith
U. Northumbria

Prototype: As few as 1.7 aerosolized infectious virions detected per min.

Efficient capture of aerosolized virus, halves infection risk.



GFP-expressing VSV-pseudotyped lentivirus nebulized using PARI TurboBoy SX and captured using PBM-HALE™ over 15 min). Condensates seeded on 10,000 HEK-293T's at 0.01 MOI. GFP expression measured at 72hrs by FACS, visualized by fluorescent microscopy. Bar = 100 µm

External validation: Efficient Aerosolised SARS-CoV-2 Particle Capture

PariBoy Classic nebuliser:


- Mean droplet diameter 3.5 μm
- 67% of mass in $< 5 \mu\text{m}$
- 5 min sampling

Particle types:

- Polystyrene beads (118 nm, -71 mV, diH₂O)
- Neutral liposomes (168 nm, -20 mV, PBS)
- Negative liposomes (188 nm, -72 mV, PBS)
- Lentiviral VLP (MLV; 193 nm, -31 mV, PBS)
- SARS-CoV-2 VLPs (155 nm, -17 V, PBS)



External validation: Efficient Aerosolised SARS-CoV-2 Particle Capture



PariBoy Classic nebuliser:

- Mean droplet diameter 3.5 μm
- 67% of mass in $< 5 \mu\text{m}$

Particle types:

- Polystyrene beads (118 nm, -71 mV, diH₂O)
- Neutral liposomes (168 nm, -20 mV, PBS)
- Negative liposomes (188 nm, -72 mV, PBS)
- Lentiviral VLP (MLV; 193 nm, -31 mV, PBS)
- SARS-CoV-2 VLPs¹ (155 nm, -17 V, PBS)



Sampling conditions:

- 5 ml sample
- 5 min nebulization
- ~1.5 ml of dry ice condensate captured

Sample analysis:

- Particle size and concentration
- Fluorescent Nanoparticle tracking analysis (Malvern; 488 nm FluoSpheres, TopFluor liposomes, YFP, or scatter)
- Unfrozen input, Frozen input, Condensate.

External validation: Efficient Aerosolised SARS-CoV-2 Particle Capture



Virtually no loss of VLP size or structure

Aggregation with highly charged molecule

- Drop in concentration
- Rise in particle size
- Beads and liposomes
- Liposome vs VLP stability?

Prototype: Clinical Pilot Update.



NO FALSE POSITIVES:

- COVID19 patients (n=12).
 - Nasal swab negative.
 - Week 2-3 of symptoms, known nasal -ve period.
 - Antivirals / hydroxychloroquine.
 - (dyspnoea) 25-30 breaths/min.
 - 5-20 min sampling.
 - In COVID19 wards.
 - Blinded analysis.
- **N=60 study ongoing**
 - Actively recruiting.
 - Must be in 0-5 days from symptom onset
 - Must be nasal positive by RT-PCR
 - N=18 completed
 - Interim data release: n=30
 - Limited attendance at week 1 of symptoms slows recruitment.

We believe we can detect COVID-19 possibly with a 1 min sample




- **Safely:** Kill the virus.
- **More reliably:** Larger sample than nasal swabs.
- **Simply:** With no skills needed: just breathe out.
- **Using mass screening:** by mass production of plastic.
- **Where patients are:** using *any* point of need testing system.
- **With current gold standard tests:**
 - e.g. Any laboratory RT-PCR instrument
 - e.g. Abbott ID NOW[®]: 5 min test or LFT/LFD
 - e.g. DeepVerge MicroTox BT: 4 sec test.

We believe we can detect COVID-19 possibly with a 1 min sample



- Confirm infectious virus load
 - By source of virus (oral, lung, nose)
 - By particle size (droplet, aerosol)
 - Optimise sampling maneuver.
- Determine the smallest sample amount needed for RT-PCR detection.
- Expand to pre/asymptomatic contacts / time course
- Deliver Emergency Use Authorisation (USA, UK).
- Produce >50,000 units.
- Supply at no profit basis under development funding.

How do I use the device?




1. Device use SOP: <https://youtu.be/h6tLt9u-rWU>
2. Lay explanation of use: https://youtu.be/TkQEj-KN_os



Dr. Pete Hotten
Chief Executive Officer
p.hotten@pulmobiomed.com



Dr. Sterghios A. Moschos FRSC FIBMS
Chief Scientific Officer
s.a.moschos@pulmobiomed.com



Scientific and Technical Team

Saqib Ali (Lead Design Engineer): Design, modification, assembly, production oversight.
Dr Theodora Mantso (Biologist): Device testing and wet biology, microbiology.

PulmoBioMed™ Ltd. is a technology spin-out of Northumbria University

Company no. 12552857

This slide deck is live and will be updated on an ASAP basis.

Acknowledgements

Moschos lab Northumbria

Dr Pep Canyelles-Pericas
Dr John Henderson
Dr Andrew Nelson
Dr William Cheung
Dr Louise Usher
Mr Kavith Shah
Mr Guilherme M.E. Silva
Mr Craig Clements
Mr Declan Gardner
Miss Zoe Hewitson

PulmoBioMed Ltd.

Dr Huw A. Edwards
Dr Pete Hotten
Mr Jonathan Brooks
Dr Theodora Mantso
Mr Saqib Ali

University of Athens

Dr Gkikas Magiorkinis
Dr Paraskevi A. Katsaounou
Mr Edison Jajah
Mr Nikolaos Athanasiou

University of Ulm

Prof Jan Münch
Dr Janis Mueller
Mr Rudiger Gross
Prof Manfred Frick

University of Oxford

Prof Aris Katzourakis

University of Crete

Prof Aristeides Tsatsakis
Prof George Sourvinos
Dr Alexandros Zafiropoulos
Dr Diamantis P. Kofteridis
Mr Taxiarchis Nikolouzakis



NATIONAL AND KAPODISTRIAN
UNIVERSITY OF ATHENS

